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1. Palmer, Robert J. (2001). "Polyamides, Plastics". *Polyamides, Plastics. Encyclopedia of Polymer Science and Technology* (4th ed.). John Wiley & Sons, Inc. doi:10.1002/0471440264.pst251.
2. [Beginners Guide to Size Exclusion Chromatography <https://www.waters.com/nextgen/global/education/primers/beginners-guide-to-size-exclusion-chromatography.html>](https://www.waters.com/nextgen/global/education/primers/beginners-guide-to-size-exclusion-chromatography.html) .
3. Wudy, Katrin & Drummer, Dietmar. (2018). Aging effects of Polyamide 12 in Selective Laser Sintering: Molecular Weight Distribution and Thermal Properties. *Additive Manufacturing* .25.10.1016/j.addma.2018.11.007.
4. Chen, An-Liu & Wei, Kuan-liang & Jeng, Ru-Jong & Lin, Jiang-Jen & Dai, Shenghong. (2010). Well-Defined Polyamide Synthesis from Diisocyanates and Diacids Involving Hindered Carbodiimide Intermediates. *Macromolecules*.44.10.1021/ma1022378.
5. [APC System Chemical Compatibility Guide. <https://www.waters.com/waters/support.htm?lid=135084097&lcid=135084096&type=USRM>](https://www.waters.com/waters/support.htm?lid=135084097&lcid=135084096&type=USRM)
6. Christian Wold, Elena Uliyanchenko (LCGC Europe), Ultra-High Size-Exclusion Separations of Engineering Plastics: Challenges and Opportunities, https://www.waters.com/webassets/cms/library/docs/lcgc_ulyanchenko_engineering_plastics_with_uhpsec.pdf <
https://www.waters.com/webassets/cms/library/docs/lcgc_ulyanchenko_engineering_plastics_with_uhpsec.pdf> .
7. Don Trinite, Waters Corp., 1360 N. Wood Dale Rd., Suite C, Wood Dale IL 60191.
8. Browne J., TA Instruments Application of Polyamides Containing Recyclates, New Castle, DE, December 2022.

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[ACQUITY UPLC RI检测器 <https://www.waters.com/134726507>](https://www.waters.com/134726507)

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© 2023 Waters Corporation. All Rights Reserved.

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